



EAGLE APPRENTICE

Recommended for: First-Time Wearers. Patients with Trouble Adapting to Other Progressives

PRIMARY PAIR

Some patients have a hard time getting used to progressive lenses. This can be especially true if they put on a pair intended for someone used to wearing a progressive lens. The Eagle Apprentice is designed to have smoother, more comfortable transitions between the far, intermediate, and near viewing zones. This makes it easier for a patient to get used to wearing a progressive. Since there's much less of a jump between these zones, they experience less swim and sway in this lens.



AVAILABILITY

MATERIALS	POWER	ADD	FITTING HEIGHTS
<ul style="list-style-type: none"> Plastic Trivex Polycarbonate 	-20.00 to +12.00	+0.75 to +4.50	14mm
<ul style="list-style-type: none"> 1.56 1.60 1.67 1.74 	CYLINDER	PRISM	16mm
	0.00 to -12.00	0.00 to 7.00	18mm
			20mm
			22mm

Exact power range availability will vary depending on prescription, material, and frame measurements. Full power range not available in all materials

COLORS

- Clear
- DayNite - Gray, Brown
- Polarized - Gray, Brown, G-15
- LifeRx - Gray, Brown
- Transitions - Gray, Brown, XTRActive, Vantage
- Photofusion - Gray, Brown
- Drivewear

Not all colors available in all materials.

FEATURES

Easy Zone Transition

Perfect for patients that have trouble adapting to progressives

Variable Decentration

Accomodate large frames and PDs that most other progressives can't handle.

Variable Inset

Optimized near viewing area tailored to each patient.

Minimized Oblique Astigmatism

Clear vision across the entire surface of the lens.

Lenticularization

Minus lenses in wrap frames are no longer a problem.

Advanced Frame Measurements

Pantoscopic Tilt

How vertically straight or tilted the frame is when worn.

Vertex Distance

Distance between eye and lens.

Wrap Angle

Degree to which the frame curves to fit the patient's face.

